
ORDERING INFORMATION			
Product Number	Description	Quantity	
849	IV, Arterial and Mid-line Catheter Securement Kit Contents: 1 securement device 1 small transparent dressing (2.375 in x 2.75 in) 1 double tape strip	50	
850	IV, Arterial and Mid-line Catheter Securement Device Contents: 1 securement device 1 double tape strip	50	
851	IV Start Kit Contents: 1 securement device 1 small transparent dressing (2.375 in x 2.75 in) 1 double tape strip 1 Prevantics® pad 1 tourniquet 2 gauze sponges (2 in x 2 in)	50	
852	CVC, PICC and Arterial Sheath Securement Device Contents: 1 securement device 1 double tape strip	50	
853	CVC, PICC and Arterial Sheath Securement Device Kit Contents: 1 securement device 1 large transparent dressing (4 in x 4.75 in) 1 double tape strip	50	
854	CVC, PICC and Arterial Sheath Securement Device Kit Contents: 1 securement device 1 large Prevahex™ CHG transparent dressing (4 in x 4.75 in) 1 double tape strip	50	
950	General Purpose Small Securement Device Contents: 1 securement device 1 double tape strip	50	
951	General Purpose Medium Securement Device Contents: 1 securement device 1 double tape strip	50	
952	General Purpose Large Securement Device Contents: 1 securement device	50	



800-343-3980

www.dalemed.com







Now You Can Get Comfortable with Security

Keeping catheters safe and secure is vital to you and your patient. A secure catheter helps prevent pistoning and accidental dislodgement while enhancing patient comfort. Dale Hold-n-Place securement dressings are designed specifically to keep the catheter positioned properly and accommodate each patient. Once the device is in place, you can be confident the line is reliably secure.

All The Features You Want

- Soft and flexible design with no hard plastic parts for improved patient comfort
- No skin prep required for application and no alcohol required for removal
- Provides superior securement for both horizontal and vertical lifting accidental line pulls
- Hypoallergenic, breathable and not made with natural rubber latex to reduce the risk of allergic reactions and skin irritation
- Protects skin from hard plastic parts
- Secures catheters and lines from the top, bottom and sides

Designed for Virtually Any Line

Seeing is Believing

We'll help you put Dale Hold-n-Place Catheter Securement to the test. Ask your Dale representative about a no-obligation, free trial offer, or call 800-343-3980.

Hold-n-Place for IV, Arterial and Mid-line Catheter Securement

- Secures IV lines and Luer hubs using the familiar chevron technique
- Maintains optimal catheter insertion angle
- Protects skin from hard plastic parts













Hold-n-Place for CVC, PICC and Arterial Sheath Securement

- Secures CVC and PICC lines using the familiar chevron technique
- Suture-free securement for protection from needlestick injuries







Hold-n-Place for General Purpose Securement

- Secures catheters and lines using repositionable hook and loop closure
- Available in small, medium and large







Dale® Nasogastric Tube Holder

The Dale Nasogastric Tube Holder is ideal for holding most nasogastric tubes, including nasoenteric, Salem Sump®, and Levin stomach tubes.

- Skin-friendly adhesive stays in place for up to 3 days, yet is easy to remove.
- Breathable material.
- Unique dual interlocking tabs spiral around the nasogastric tube for superior tube security.
- Non-adherent plastic tips on both tabs make tube removal safer and easier, especially when wearing gloves.
- Single patient use only.
- No scissors required.
- May also be used to secure most nasal catheters used for manometry and pH impedance testing





ORDERING INFORMATION		
Product Number	Size	Quantity
160	Large / Adult	50 per box
161	Small / Young Adult / Pediatric / Geriatric	50 per box

IV Line Solutions

Dale® IV-ARMOR®

Dale's IV-ARMOR minimizes downstream occlusions while enabling patient movement. The flexible protective overlay is the ideal way to maximize IV patency while protecting the site from patient tampering.

- Reduces the introduction of infection associated with IV reinsertion.
- Reduces downstream occlusion alarms.
- Reduces risk of catheter kinking in a bend such as a wrist or elbow.
- Reduces risk of catheter movement in an area of bony prominence, such as back of hand or forearm.
- Skin friendly adhesive.









ORDERING INFORMATION		
Product Number Size Quantity		
675	One size fits most	10 per box

Dale® Bendable ArmBoard

The Dale Bendable ArmBoard allows the patient's hand and wrist to maintain their natural position while securely and comfortably preventing catheter movement.

- Bendable inner core can be custom shaped to desired position - Enables clinician to hyperextend patient's wrist for pronounced radial artery exposure, trouble-free insertion and line integrity.
- Adjustable hook & loop closures provide comfortable, snug fit.
- Soft, plush lining enhances patient comfort.
- Four convenient sizes.
- MR Conditional.







ORDERING INFORMATION		
Product Number	Size	Quantity
650	Large, 23cm x 9cm (9"x3.5")	10 per box
651	Medium, 13.7cm x 6.4cm (5.38"x2.5")	10 per box
652	Small, 11.4cm x 4.4cm (4.5"x1.75")	10 per box
653	X-Small, 10.8cm x 2.5cm (4.25"x1")	10 per box

Respiratory

Dale® Stabilock Endotracheal Tube Holder

The Dale Stabilock Endotracheal Tube Holder provides a fast and easy way to secure endotracheal tubes and help prevent accidental extubation.

- Cushioned neck band is soft and comfortable on the face and neck.
- Skin-friendly adhesive can be left in place up to 3 days.
- Tubing channel is made of lightweight, flexible material to securely hold the endotracheal tube.
- Hook & loop closures provide a secure fit.
- Easy repositioning reduces risk of skin and lip ulcerations.
- Easy access for oral care.





ORDERING INFORMATION		
Product Number	Size	Quantity
270	Fits endotracheal tube sizes 7.0mm-10.0mm	10 per box
273	Adhesive base	10 per box

Dale® Tracheostomy Tube Holders

The Dale Tracheostomy Tube Holder is designed to provide secure positioning and minimize movement of the tracheostomy tube. Unnecessary movement of the trach tube can result in accidental dislodgement or displacement of the tube, tracheal fistula, tracheal stenosis, or airway granuloma.

- Eliminates the frustrations associated with twill ties and other holders while minimizing secondary complications.
- Moisture wicking neck band lining reduces the risk of skin breakdown.
- Fastener tabs quickly and securely attach to the tracheostomy tube holder.









ORDERING INFORMATION			
Product Number	Product Number Size		
240	240 Blue® 2.54cm(1") wide neckband, fits up to 49.5cm(19.5") neck	10 per box	
241	241 PediPrints® 2.54cm(1")wide neckband, fits up to 46cm(18") neck	10 per box	
242	242 PediPrints® 1.91cm(0.75") wide neckband, fits up to 23cm(9") neck	10 per box	
244	244 Blue $^{\! 8}$ for bariatric patients 3.81cm(1.5") wide neckband, fits up to 71cm(28") neck	10 per box	
200	Dale® One Piece Holder 2.54cm(1") wide neckband, fits 24.1cm-45.5cm(9.5" - 18") neck	30 per box	

Medical/Surgical

Dale® Post-Surgical Bra

The Dale Post-Surgical Bra provides the support and compression patients need immediately following breast surgery. Soft, breathable fabric stretches to accommodate edema without binding or irritating sensitive tissue.

- Front hook & loop closure facilitates inspection of the surgical site and dressing changes.
- Ideal for use during radiation treatment.
- Can be worn alone or can be used to hold surgical dressings in place without tape.
- Indicated for use after a wide range of surgical procedures: Reduction, reconstruction, lift, augmentation, mastectomy, lumpectomy, biopsy, and fine needle aspiration.






ORDERING INFORMATION		
Product Number	Size	Quantity
702	Medium, fits B-D 86cm-92cm(34"-36")	1 per box
703	Large, fits B-D 92cm-97cm(36"-38")	1 per box
704	X-Large, fits B-D 97cm-112cm(38"-44")	1 per box
705	XX-Large, fits C-E 117cm-137cm(46"-54")	1 per box

Dale® Abdominal Binder with EasyGrip™ Strip

Dale Abdominal Binders are the only binders clinically shown to expedite the return of pulmonary function after open abdominal surgery while encouraging patients to ambulate sooner. Available in sizes up to 94" long and in widths up to 15".

- Exclusive stretch material provides full, all-around compression.
- Hook & loop closure secures at any point along the binder circumference.
- Promotes mobility and lessens the chance of sedentary complications.
- Available in bariatric sizes.
- EasyGrip[™] Strip holds up to four 100cc drainage bulbs for easier monitoring and changing.







ORDERING INFORMATION		
Product Number	Size	Quantity
525	6" wide, stretches to fit 72cm-127cm(28"-50")	1 per box
410	9" wide, stretches to fit 76cm-114cm(30"-45")	1 per box
411	9" wide, stretches to fit 117cm-157cm(46"-62")	1 per box
418	9" wide, stretches to fit 152cm-191cm(60"-75")	1 per box
810	12" wide, stretches to fit 76cm-114cm(30"-45")	1 per box
811	12" wide, stretches to fit 117cm-157cm(46"-62")	1 per box
818	12" wide, stretches to fit 152cm-191cm(60"-75")	1 per box
820	12" wide, stretches to fit 183cm-239cm(72"-94")	1 per box
920	15" wide, stretches to fit 183cm-239cm(72"-94")	1 per bag

Medical/Surgical

Dale® Nasal Dressing Holder

The Dale Nasal Dressing Holder replaces adhesive tape to hold nasal dressings in place following rhinoplasty, septoplasty, sinus surgery, facial trauma and nose bleeds. Clear plastic dressing holder is secured by ear loops. Hook and loop fastener tabs can be adjusted to fit any patient comfortably.

- Soft moisture repellent material adds to comfort.
- Plastic center section holds dressings in place without absorbing secretions.
- Easy to remove and replace, making dressing changes simple and painless.



ORDERING INFORMATION		
Product Number	Size	Quantity
600	One size fits most	10 per box

Dale® Transducer Holder

The Dale Transducer Holders offers a level and secure platform for up to three transducers. Stability ensures correct waveform readings and corresponding numerical display.

- Platform conforms to patient's arm.
- Soft, comfortable material with adjustable hook and loop closures.
- Three holding tabs reduce the time associated with manual re-leveling and taping individual transducers to the patient's chest, arm, or mounting on the bedside pole.
- Increases organization and security of hemodynamic catheters during transportation.




ORDERING INFORMATION		
Product Number Size Quantity		
620	One size fits most	10 per box

Urology

Dale® Hold-n-Place® Foley Catheter Holders

Dale Hold-n-Place Foley Catheter Holders stabilize the indwelling catheter, reducing the likelihood of traumatic removal, without restricting patient's movement or irritating the skin. All Dale Hold-n-Place Foley Catheter Holders feature unique locking systems that securely anchor the catheter in place, reducing the risk of urethral erosion, bladder spasms and trauma.

- The leg band, made entirely of soft stretch material, distributes compression evenly around the thigh to minimize circumferential pressure. Ideal for long-term catheterization.
- The adhesive patch, with a skin-friendly adhesive base, is ideal for short-term catheterized patients.
- The locking devices speed up application and simplify adjustment.
- The holders can be used with equal effectiveness on all Foley catheters including Teflon, silver and silicone coated catheters.
- Leg band/waist band ideal for bariatric patients.



Hold-n-Place Adhesive Patch







ORDERING INFORMATION		
Product Number	Size	Quantity
316	Leg Band, fits up to 76cm(30")	10 per box
330	Leg Band/Waist Band Fits up to 142cm(56")	10 per box
150	Adhesive Patch, One size fits most	50 per box
490	Urodynamic Holder, One size fits most	50 per box