




RED CELL DOCKABLE FILTRATION 0459

Ref. KFS0001XB



SPECIFICATIONS

Technical specifications

Product code	Filter	Packaging	
KFS0001XB	LEUCOBED LCG3	1 device / peelable sachet	12 devices / box

Included items

- 1 Spike
- Air-vents
- Drip chamber
- Pre-filter (200 µm)
- Rolled tubing
- Roller clamp

Sterilization

β Irradiation

Instructions for use

NOTEURO11

Indications

The MacoPharma LeucoBed LCG3 High Efficiency Leucocyte Removal Filter is a bedside filter intended for the removal of leucocytes, platelets and microaggregates from a single unit of red cells / whole blood.

Sterilized by β Irradiation.

Sterile and non-pyrogenic fluid pathway.

Precautions

Follow instructions carefully

Follow normal aseptic techniques

For single use : do not re-sterilize or re-use

Do not use if packaging is damaged or end protectors are lost or displaced

For a better filtration and a better storage of filtered red cells or buffy-coat removed red cells, it is recommended that the blood used is not older than 15 days.

Use

Leucocytes Removal Filter for Red Cell Concentrate or Whole blood

Label

English

Shelf Life

3 years