




# CORD BLOOD

## COLLECTION BAG BY MACO PHARMA

Maco Pharma's cord blood range provides products to collect and to cryopreserve stem cells of the umbilical cord blood.

As one of the world's first companies to provide cord blood banks with dedicated collection kits, Maco Pharma has been a leader in the field for more than fifteen years. Based on the company transfusion deep knowledge, collection bags with one or two needles were designed to collect cord blood from both caesarean section deliveries and natural births.



## CORD BLOOD COLLECTION




Maco Pharma collection bags are specifically designed for the collection of human umbilical cord blood. The collection bags are sterile and suitable for vaginal births and caesarean section deliveries in operating theaters.

### Cord blood Collection Bags

References	Nominal volume (mL)	Collection volume (mL)	Dimensions (cm) Tubing length	Anticoagulant solution	Connections	Collection lines	Needle size	Sterilisation	Packaging	Shelf life
MSC1200PU	250	200 mL when the additional 8 mL of CPD* have been transferred from the rinsing bag to the collection bag	150	21 mL CPD* in the collection bag + 8 mL CPD in the rinsing bag	1 outlet port 1 injection site	1	12 G	Steam	Double peelable overwrap	20/ box
MSC1210PU		2								
MSC1230PU		150 mL		21 mL CPD* in the collection bag	1					
MSC1206DU	300	200 mL when the additional 8 mL of CPD* have been transferred from the rinsing bag to the collection bag	152	21 mL CPD* in the collection bag + 8 mL CPD in the rinsing bag	2 outlet ports 1 injection site	1				
MSC1208DU		2								

\*Citrate Phosphate Dextrose



## MACO BIOTECH FREEZING - EVA BAGS

CE 0459



## EVA Bag for cell freezing - Mini 3D Bag

Reference	Volume (mL)	Dimensions (mm)			Connections	Sterilisation	Packaging		Conditions of use	Shelf life
		Length	Width*	Height						
GSR0501AU	25 mL (A: 20 mL + B: 5 mL)	86	65.5	8	1 spike 1 male luer lock 2 needleless access connectors 2 outlet ports	Beta Irradiation	Single peelable overwrap	25/box	Temperature of use: down to -196°C Temperature of thawing: +37°C	5 years

\*With the outlet ports

## Protective container for freezing - Mini 3D Overwrap

Reference	Dimensions (mm)			Compatible with	Sterilisation	Packaging		Conditions of use	Shelf life
	Length	Width	Height						
GSR0010AU	92	101.5	9.5	GSR0501AU	Beta Irradiation	Single peelable overwrap	75/box	Temperature of use: down to -196°C Temperature of thawing: +37°C	2 years



GSR0501AU



GSR0010AU