


Filtration of Red Cell Concentrate









Leucoflex LCR5
Integrated to a
Quadruple bag
CPD/SAGM, 450ml

Ref. LQE 6285PQ

Specifications

- . **Filter** : Leucoflex LCR5
- . **Bags** : 450ml PVC
2 Transfert bags made in TOTM
- . **Solutions** : 63ml CPD, 100ml SAGM
- . **Supplementary devices** : Vacuam,
Secuam
- . **Labels** : four-language label

Use




Packaging

1 bag / peelable sachet - 12 bags / box



Whole Blood Filtration



Leucodepleted Red Cell Concentrate
Leucodepleted Plasma

Numbered Tubing

Standard bag

Platelet Storage bag (TOTM)





0459

Ref. DQE6281LQ

Specifications

- . **Filter** : LEUCOFLEX LST1
- . **Bags** : 4 PVC
- . **Solutions** : 63 mL CPD / 100 mL SAGM
- . **Collection** : 450 mL +/- 10 %
- . **Included items** : Bactivam, Vacuvam, Secuvam
- . **Labels** : English, French
- . **Sterilisation** : Steam
- . **Shelf life** : 2 years
- . **Packaging** : 1 pack/peelable sachet - 12 packs/box

Use



■ Whole Blood
 ■ RCC
 ■ PLASMA
 ■ MB PLASMA
 ■ PRP
 ■ Platelets
 ■ SAG-M PAGGSM
 ■ NaCl
 ■ Buffy Coat
 ■ Serum
 ■ Supernatant
 ■ Cryoprecipitate

Certified ISO 9001
Certified ISO 13485

MacoPharma
F - 59420 MOUVAUX

DQE6281LQ_A_A 07/09/2006